



COBA Care Disinfectant Solution

CaterQAC, a quaternary ammonium based disinfectant, is highly recommended for use against Covid-19

- Ideal disinfectant solution for use with disinfectant mats, to sanitise footwear.
- Application for disinfectant mats:
- Add recommended amount to the mat basin. Step onto the mat, covering base of footwear in dilution. Gently scrape and sanitise by wiping soles on the mat surface. Renew the solution when visibly dirty.
- Use with COBA Disinfectant Mat or Fingertip Disinfectant Mat
- For use on all surfaces, in all environments.
- Food approved – SANS 1828.
- PH in dilution: 10.5 +- 0.5.
- Active ingredient: Di-n-decyl dimethyl ammonium chloride.
- Avoid contact with skin and eyes (see product label for further precautions).
- Strength designation: 4.6%.
- Dilution Rate: 300ml to 10L water



Parts

Part Number	Size	Colour	Weight (kg)
CATERQAC-20LT	20 L		
CATERQAC-5LT	5 L		

Technical Specifications

Max. Operating Temperature	Store below 30°C.
Typical Applications	Ideal disinfectant solution for use with disinfectant mats
Other Product Performance Information	PH in dilution: 10.5 +- 0.5.
Antibacterial effectiveness (GB/T 31402-2015)	Strength designation: 4.6%.

Additional Information



Get in touch | To order FREE samples or arrange a site visit please call +27 (0) 11 452 7961 or email sales@cobaafrika.com

www.matting.co.za